INSTRUCTIONS FOR USE














Ambu

Ambu® aScope™ Gastro

For use by trained clinicians/physicians only.
For use in medical facility.
For use with Ambu® aBox™ 2.

QUICK GUIDE

Pat. Pending

2797

Ambu is a registered trademark and aScope and aBox are trademarks of Ambu A/S.

Contents

Page

English (Instructions for use).....	4-20
Български (Инструкция за употреба)	21-40
Čeština (Návod k použití)	41-58
Dansk (Brugervejledning)	59-75
Deutsch (Bedienungsanleitung)	76-94
Ελληνικά (Οδηγίες χρήσης).....	95-114
Español (Instrucciones de uso).....	115-134
Eesti (Kasutusjuhend)	135-151
Suomi (Käyttöohje)	152-169
Français (Mode d'emploi).....	170-187
Hrvatski (Upute za upotrebu).....	188-205
Magyar (Használati útmutató).....	206-223
Italiano (Istruzioni per l'uso).....	224-241
日本語（使用説明書）	242-257
Lietuviškai (Naudojimo instrukcija).....	258-275
Latviski (Lietošanas instrukcija)	276-293
Nederlands (Gebruiksaanwijzing).....	294-312
Norsk (Bruksanvisning).....	313-329
Polski (Instrukcje użytkowania).....	330-347
Português (Instruções de utilização)	348-365
Română (Instrucțiuni de utilizare)	366-383
Русский (Инструкция по применению).....	384-404
Slovenčina (Návod na použitie)	405-422
Slovenščina (Navodila za uporabo)	423-440
Svenska (Bruksanvisning)	441-457
Türkçe (Kullanma talimatları)	458-475
中文 (使用说明)	476-491

1. Important Information – Read Before Use!

Read the *Instructions for Use* (IFU) carefully before using the aScope™ Gastro. These instructions describe the function, set-up and precautions related to operating the aScope Gastro. Please be aware that these instructions do not describe clinical procedures. Prior to use of the aScope Gastro, it is important for operators to have received sufficient training in clinical endoscopic techniques and to be familiar with the intended use, warnings, cautions, indications and contraindications mentioned in these instructions. There is no warranty for the aScope Gastro. In this document aScope Gastro refers to instructions which apply to the endoscope only, while *system* often refers to information relevant for the aScope Gastro and the compatible aBox™ 2 displaying unit and accessories. The *Instructions for Use* may be updated without further notice. Copies of the current version are available upon request.

1.1. Intended Use/Indications For Use

The aScope Gastro is a sterile, single-use, flexible gastroscope intended to be used for endoscopic access to and examination of the upper gastrointestinal anatomy. The aScope Gastro is intended to provide visualization via a compatible Ambu displaying unit and to be used with endotherapy accessories and other ancillary equipment.

1.2. Intended Patient Population

The aScope Gastro is intended to be used in adults; this means patients with an age of 18 years or above. The aScope Gastro is used in patients with indications in the upper gastrointestinal anatomy requiring visualization and/or examination with flexible gastroscopy and use of endoscopic accessories and/or equipment.

1.3. Contraindications

No known contraindications.

1.4. Clinical Benefits

The aScope Gastro when used with the compatible aBox 2 displaying unit enables visualization, examination and endoscopic intervention of key anatomical structures in the upper GI tract, particularly esophagus, gastroesophageal junction, stomach, pylorus, duodenal bulb and descending duodenum.

High-definition imaging technology will enable endoscopists to view finer mucosal and vascular details.

The risk of patient cross-contamination is eliminated, as compared to reusable endoscopes, as the aScope Gastro is a sterile single-use medical device. Also, the risk of anaphylactic shock associated with exposure to high level disinfectants (as is possible when using reprocessed scopes) is eliminated.

1.5. Warnings and Cautions

Warnings

1. For single-use only. Do not reuse, reprocess or resterilize as these processes may leave harmful residues or cause malfunction of the aScope Gastro. Reuse of the aScope Gastro may cause cross-contamination potentially leading to infections.
2. Confirm that the outlet of the insufflation/rinsing valve is not blocked. If the opening is blocked, gas is fed continuously and patient pain, bleeding, perforation and/or gas embolism can result.
3. Prior to use always perform an inspection and functionality check according to sections 4.1 and 4.4. Do not use the device if the aScope Gastro or its packaging is damaged in any way or if the functionality check fails, as this can lead to patient injury or infection.
4. Patient leakage currents may be additive, when using energised endoscopic accessories. Do not use energised endoscopic accessories which are not classified as "type CF" or "type BF" applied part according to IEC 60601, as that could lead to too high patient leakage current.
5. Do not perform procedures with High Frequency endoscopic accessories if flammable or explosive gases are present in the gastrointestinal tract as this may result in serious injury to the patient.
6. Always observe the live endoscopic image when inserting, withdrawing or operating the aScope Gastro. Failure to do so may result in patient injury, bleeding and/or perforation.
7. Ensure that the insufflator is not connected to the auxiliary water inlet as this may cause overinsufflation which can result in patient pain, bleeding, perforation and/or gas embolism.
8. The distal end of the aScope Gastro may get warm due to heating from the LEDs. Avoid long periods of contact between the distal end of the aScope Gastro and the mucosal membrane as sustained contact may cause mucosal injury.
9. Do not insert or withdraw the aScope Gastro if endoscopic accessory is protruding from the distal end of the working channel as this may result in injury to the patient.
10. If the biopsy cap is left uncapped or if the biopsy valve is damaged it can reduce the efficacy of the aScope Gastro's suction system, and may leak or spray patient debris or fluids, posing a risk of infection. When the valve is uncapped, place a piece of sterile gauze over it to prevent leakage.
11. Always use gauze to pull the accessory device through the biopsy valve as patient debris or fluids may leak or spray, posing a risk of infection.
12. During the procedure always wear personal protective equipment (PPE) to protect against contact with potentially infectious material. Failure to do so may cause contamination potentially leading to infections.
13. Using HF (High Frequency) tools/electrosurgical equipment with aScope Gastro may disturb the image on the displaying unit which may lead to patient injury. To reduce disturbance try alternative settings on HF generator with lower Peak Voltage.
14. Portable RF (Radio Frequency) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the aScope Gastro and the displaying unit, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result which could lead to patient injury.

Cautions

1. Only use the aScope Gastro with medical electrical equipment that complies with IEC 60601-1, and any applicable collateral and particular standards. Failure to do so may lead to equipment damage.
2. Prior to using any high frequency accessory device, check the compatibility with the aScope Gastro. Always follow the instructions for use for the third party device. Failure to do so may lead to equipment damage.
3. Do not activate energized endoscopic accessories before the distal end of the endoscopic accessory is in the field of view and is extended at an appropriate distance from the distal end of the endoscope as this may result in aScope Gastro damage.
4. Do not apply oil based lubrication in the working channel as this may increase friction when inserting instruments.
5. Do not coil the insertion tube or umbilical cord to a diameter of less than 20 cm as this may damage the aScope Gastro.
6. Excessive force – Do not drop, bump, bend, torque or pull any portion of the aScope Gastro with excessive force as the aScope Gastro may get damaged leading to failure in functionality.
7. Do not use excessive force to advance the accessory device through the working channel. Doing so may cause damage to the working channel of the aScope Gastro.

1.6. Potential Adverse Events

Possible complications include (not exhaustive):

- Air embolism
- Gagging
- Gastric-to-pulmonary aspiration
- Mucosal laceration
- Mucosal bleeding
- Perforation
- Peritonitis

1.7. General Notes


If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

2. Endoscope Description

aScope Gastro must be connected to the aBox 2 displaying unit. Please refer to the aBox 2 *Instructions for Use* for more set-up details.

2.1. Description of Components

aScope Gastro

Ambu® aScope™ Gastro Sterile and Single-Use Gastroscope	Part Number	Distal End Outer Diameter	Working Channel Inner Diameter
	483001000	9.9 mm	2.8 mm

2.2. Description and Function

The aScope Gastro is a sterile and single-use gastroscope for use within the upper gastrointestinal tract. The aScope Gastro is inserted into the patient's upper gastrointestinal tract through the mouth. The aScope Gastro is powered by connection to the aBox 2 displaying unit. The aScope Gastro can be used with endotherapy accessories and ancillary equipment for endoscopic procedures. The working channel allows for the passage of endoscopic accessories, instillation of fluids, and suction of fluids.

The Water Jet function allows the instillation of fluids. The insufflation function allows the instillation of CO₂ to expand the gastrointestinal lumen. Insufflation has multiple purposes hence specification may not be necessary.

The aScope Gastro is intended for left-handed use.

The optical module in the distal tip consists of a camera housing which contains a camera and LED light sources.

The user can angulate the distal tip in multiple planes for visualization of the upper gastrointestinal tract by turning the control wheels to activate the bending section (see figure 1). The bending section can bend up to 210° enabling a retroflexion to visualize the fundus and esophageal sphincter.



Figure 1

No.	Part	Function
	Control Section	The user holds the aScope Gastro at the control section with the left hand. The control wheels and remote switches may be manipulated with the left and/or right hand.

No. on Fig.1	Part	Function
1	Remote Switches/ Programmable Buttons	The user activates functions on the aBox 2. The functions of the remote switch buttons are pre-configured from factory and can be re-configured according to the user's preference.
2	Up-Down Control Wheel	The Up/Down control wheel manipulates the bending section of the scope. When this wheel is turned in the "U" direction, the bending section moves UP; when the wheel is turned in the "D" direction, the bending section moves DOWN.
3	Up/Down Angulation Lock	Turning this lock in the "F" direction frees angulation. Turning the lock in the opposite direction locks the bending section at any desired position along the up/down axis.
4	Suction Valve	The removable suction valve activates suction when pressed down. The valve is pressed to activate suction to remove any fluids, debris or gas from the patient.
5	Insufflation & Rinsing Valve	The insufflation/rinsing valve controls insufflation and lens rinsing. Placing a finger on the opening of the valve activates insufflation. The valve activates lens washing when pressed down.
6	Right-Left Control Wheel	The Right/Left control wheel manipulates the bending section of the scope. When this wheel is turned in the "R" direction, the bending section moves RIGHT; when the wheel is turned in the "L" direction, the bending section moves LEFT.
7	Right/Left Angulation Lock	Turning this lock in the "F" direction frees angulation. Turning the lock in the opposite direction locks the bending section at any desired position along the right/left axis.
8	Biopsy Valve	Biopsy access valve. Seals the working channel.




9	Working Channel	The working channel functions as: <ul style="list-style-type: none"> • Suction channel. • Channel for the insertion or connection of endoscopic accessories. • Fluid feed channel (from a syringe via the biopsy valve).
10	Insertion Tube	The flexible insertion tube is inserted into the patient's upper gastrointestinal tract.
11	Bending Section	The bending section is the maneuverable part of the aScope Gastro, that can be controlled by the control wheels and angulation locks.
12	Distal Tip	The distal tip holds the camera, the light source (two LEDs), the working channel exit, the insufflation & rinsing nozzle, and the Water Jet outlet.
13	Insufflation & Rinsing Nozzle	Nozzle for lens rinsing & insufflation.
14	Water Jet Outlet	Water jet is used for endoscopic irrigation of the patient's upper gastrointestinal tract.
15	Working Channel Outlet	The working channel outlet allows: <ul style="list-style-type: none"> • Aspiration of fluids. • Introduction of endoscopic accessories. • Introduction of fluids.
16	Camera	Enables visualization of the upper gastrointestinal tract.
17	LEDs	Enables illumination of the upper gastrointestinal tract.
18	Suction Connector	Connects the aScope Gastro to the suction tube of the suction pump.
19	Water Jet Connector	Connects the aScope Gastro to the irrigation tube of the irrigation pump. The Water Jet connector has an integrated one-way valve to reduce the risk of backflow.
20	aScope Gastro Connector	Connects the aScope Gastro to the output socket of the aBox 2. Ancillary equipment for suction, insufflation, lens rinsing and irrigation can be attached to the aScope Gastro connector.
21	Insufflation & Rinsing Connector	Connects the aScope Gastro to the sterile water bottle to enable insufflation/lens rinsing.
22	Disconnection Button	Press the button when disconnecting the aScope Gastro from the aBox 2.
23	Spare Suction Valve	Can be used to replace the existing suction valve in case of blockage.

2.3. Product Compatibility

The aScope Gastro can be used in conjunction with:

- The Ambu aBox 2.
- Insufflator for endoscopic gastrointestinal procedures: source of medical grade CO₂.
- Standard disposable insufflation/rinsing fluid management tubing sets compatible with Olympus.
- Vacuum source to provide a vacuum.
- Standard suction tubes.
- Gastrointestinal endoscopic instruments specified to be compatible with a working channel size (ID) of 2.8 mm or less (There is no guarantee that instruments selected using only this minimum working channel size will be compatible with the aScope Gastro).
- Medical-grade water-based lubricants, iodine based contrast agents, lipiodol, haemostatic agents, lifting agents, antifoaming agents, tattoo for permanent staining, and dyes for vital staining.
- Sterile water.
- High frequency electrosurgical equipment fulfilling EN 60601-2-2. To keep high frequency leakage currents within allowed limits, the maximum sinus peak voltage level of the electrosurgical unit shall not exceed 5.0 kVp.
- Irrigation pump for endoscopic gastrointestinal procedures with luer connector.

3. Explanation of Symbols Used

Symbol	Description	Symbol	Description
	Working length of the insertion tube.	Rx only	US Federal Law restricts this device to sale by or on the order of a physician.
	Maximum insertion portion width (Maximum outer diameter).		Atmospheric pressure limitation: Between 80 and 106 kPa in operating environment.
	Minimum working channel width (Minimum inner diameter).		Humidity limitation: Relative humidity between 30 and 85 % in operating environment.
	Country of manufacturer: Made in Malaysia		Temperature limit.
	Field of view.		Medical Device.
	Warning.		Packaging level ensuring sterility.
	Do not use if package is damaged.		Global Trade Item Number.
	IFU symbol.		UL Recognized Component Mark for Canada and the United States.

The complete list of symbol descriptions can be found on ambu.com.

4. Use of the aScope Gastro

The number in grey circles refer to the Quick Guide on page 2. Before each procedure, prepare and inspect a new aScope Gastro as instructed below. Inspect other equipment to be used with this aScope Gastro as instructed in their respective instruction manuals. Should any irregularity be observed after inspection, follow the instructions as described in section 6, "Trouble Shooting". If the aScope Gastro malfunctions, do not use it. Contact your Ambu sales representative for further assistance.

4.1. Inspection of the aScope Gastro 1

- Check that the pouch seal is intact and discard the aScope Gastro if the sterile seal has been damaged. **1a**
- Carefully peel off the peel pouch packaging of the aScope Gastro and remove the protective elements from the handle and the distal end. **1b**
- Carefully run your hand back and forth over the entire length of the insertion tube, including the bending section and distal end, of the aScope Gastro to make sure that there are no impurities or damage on the product such as rough surfaces, sharp edges or protrusions which may harm the patient. Make sure to use aseptic technique when performing the above. Otherwise the sterility of the product will be compromised. **1c**
- Inspect the distal end of the aScope Gastro's insertion tube for scratches, cracks or other irregularities. **1c**
- Confirm that the top hole of the Insufflation/Rinsing valve is not blocked. **1c**
- Turn the Up-Down and Right-Left control wheels in each direction until it stops and return to neutral position. Confirm that the bending section functions smoothly and correctly, that maximum angulation is achieved and returned to the neutral position. **1d**
- Confirm that the angulation locks are functional by locking and releasing the angulation locks according to section 2.2. Turn the angulation wheels fully in all directions, lock the angulation in a fully angulated position and confirm that the bending section can be stabilized. Release the angulation locks and confirm that the bending section straightens out.
- Using a syringe, flush sterile water into the working channel. Ensure that there are no leaks and that the water is emitted from the distal end. **1e**
- If applicable, confirm compatibility with applicable accessory devices.
- A spare suction valve is available and placed securely on the mounting card together with the aScope Gastro.
- A new aScope Gastro should be readily available so the procedure can be continued in case a malfunction should occur.

4.2. Preparations for Use

Prepare and inspect aBox 2, water container, suction and irrigation pumps, and CO₂ insufflator as described in their respective instruction manuals.

- Power up the aBox 2 **2**. Carefully align the arrows on the connector of the aScope Gastro cable with the port of the aBox 2 to prevent damage to the connectors. **3**
- Connect the aScope Gastro to the aBox 2 by plugging the aScope Gastro connector into the corresponding aBox 2 female connector.
- Check that the aScope Gastro has been firmly locked to the aBox 2.
- When using the aScope Gastro it is recommended to use a mouthpiece to prevent the patient from accidentally biting the insertion tube.

4.3. Attaching Ancillary Equipment

The aScope Gastro is designed to work with most commonly available medical suction and fluid management systems.

Irrespective of the chosen fluid management system, overflow protection must be a feature of the suction canister assembly utilized in order to prevent fluids from entering the system.

This feature is commonly referred to as "self-sealing" feature or a "shut-off-filter", or similar mechanism. Please note that a new suction canister and connection should be used for each new procedure.

The aScope Gastro does not itself produce negative pressure and therefore an external vacuum source (e.g. wall suction or medical grade suction pump) will be required to operate the system. Standard suction tubings with a nominal diameter should be sufficient, as long as they simply and securely fit over the standard size suction connector of the aScope Gastro.

Nomenclature used in this section follows the established convention that each canister utilized will have multiple connection ports labeled either "To Vacuum" or "To Patient". However, it is the responsibility of the user to follow all third party manufacturer instructions and guidance applicable to the endoscopic fluid management system chosen for use with the aScope Gastro.

When using third party devices with the aScope Gastro always consult the instructions for use accompanying the third party device.

Connection to the Insufflation/Rinsing System 4

The aScope Gastro is designed to work with a source of medical grade CO₂ for insufflation. Connect the aScope Gastro using a standard disposable insufflation/rinsing fluid management tubing set.

To perform patient examinations or procedures, all fluid containers (e.g. sterile water bottle) must be properly and securely placed in order to prevent spillage, thereby maintaining a safe working environment. Place the containers in the designated locations and connect them according to the instructions in this section. Please note that a new suction canister/sterile water bottle must be used for each new procedure.

- If ancillary equipment is ON turn OFF.
- Connect the fluid container's connection adapter to the aScope Gastro.
- Confirm that the fluid container's connection adapter fits properly and that it cannot be rotated.
- Turn the ancillary equipment back ON.

Connection to the Suction System 5

Regardless of the vacuum source chosen, the aScope Gastro will require the source to provide a vacuum for the aScope Gastro to operate normally. Failure to provide the minimum vacuum requirements could result in a decreased capacity in removing patient waste and/or irrigation fluid.

- If ancillary equipment is ON turn OFF.
- Fit the end of the suction tubing securely over the suction connector located on the aScope Gastro connector.
- Connect the other end of the suction tubing from the suction connector of the aScope Gastro to the external suction system (e.g. wall suction or medical suction pump).
- Turn the ancillary equipment back ON.

Connection to the Auxiliary Water Jet System 5

- The aScope Gastro has a Auxiliary Water Jet connector with an integrated one-way valve to reduce the risk of backflow.
- If ancillary equipment is ON turn OFF.
- Just connect the **Auxiliary water** tubing to the auxiliary water inlet located on the aScope Gastro connector.
- Confirm that the fluid container's connection adapter fits properly.
- Turn the ancillary equipment back ON.

4.4. Inspection of the Endoscopic System

Checking the Working Channel 6

- Confirm that the biopsy valve is attached to the working channel port.
- Endoscopic instruments labeled for use with a working channel size (ID) of 2.8 mm or less are compatible.
- There is no guarantee that endoscopic instruments selected solely using this working channel size will be compatible with the aScope Gastro.
- Compatibility of selected endoscopic instruments should be tested prior to procedure.

Inspection of the Image 7

- Verify that a live video image and correct orientation appears on the monitor by pointing the distal end of the aScope Gastro towards an object, e.g. the palm of your hand.
- Adjust the image preferences on the aBox 2 if necessary (please refer to the aBox 2 *Instructions for Use*).
- If the image is impaired and/or unclear, wipe the lens at the distal end using a sterile cloth.
- The images must not be used as an independent diagnostic of any pathology. Physicians must interpret and qualify any findings by other means and in the light of the patient's clinical characteristics.

Checking the Remote Switches

- All remote switches should be checked to work normally even if they are not expected to be used.
- Press every remote switch and confirm that the specified functions work as expected.

Checking the Suction, Rinsing and Insufflation Functionality

- Check that the suction and insufflation/rinsing valves work as expected by pressing the suction and insufflation/rinsing valves.
- Cover the opening of the insufflation/rinsing valve and confirm that the insufflation function works properly.

Checking the Auxiliary Water Jet Functionality

- Check that the auxiliary water jet system by activating the auxiliary water pump and confirm that irrigation function works properly.

4.5. Operating the aScope Gastro

Insertion of the aScope Gastro 8

- Insert a suitable mouthpiece and place it between the patient's teeth or gums.
- If necessary, apply a medical-grade lubricant to the insertion section.
- Insert the distal end of the aScope Gastro through the opening of the mouthpiece, then from the mouth to the pharynx while viewing the endoscopic image. Do not insert beyond the proximal end maximum length mark.

Holding and Maneuvering the aScope Gastro

- The control section of the aScope Gastro is designed to be held in the operator's hand.
- The insufflation/rinsing and suction valves can be operated using a finger.
- The UP/DOWN angulation control wheel can be operated using the fingers.
- The operator's hand is free to manipulate the insertion section and the RIGHT/LEFT angulation control wheel.

Angulation of the Distal End

- Operate the angulation control wheels as necessary to guide the distal end for insertion and observation.
- The aScope Gastro's angulation locks are used to hold the angulated distal end in position.

Insufflation/Rinsing

- Cover the insufflation/rinsing valve's hole to feed CO₂ from the insufflation/rinsing valves at the distal end.
- Press the insufflation/rinsing valve to feed sterile water onto the objective lens.
- Always dispose of insufflation/rinsing tube and corresponding water bottle after each endoscopic procedure.

Instillation of Fluids

- Fluids can be injected through the working channel by inserting a syringe into the working channel port of the aScope Gastro. Insert the syringe completely into the working channel port and press the plunger to inject fluid.
- Make sure you do not apply suction during this process, as this will direct the injected fluids into the suction collection system. To ensure that all fluids have left the channel, flush the channel with approx. 5 seconds of CO₂.

Auxiliary Water Jet System

- Activate the auxiliary Water Jet system in order to apply irrigation.
- A delay in irrigation may be experienced if the auxiliary Water Jet system has not been pre-filled during the pre-procedural preparation.
- Always dispose of Water Jet tubing and corresponding water bottle after each endoscopic procedure.

Suction

- Press the suction valve to aspirate excess fluids or other debris obscuring the endoscopic image.
- For optimal suction capability it is recommended to remove accessories entirely during suction.
- Should the suction valve on the aScope Gastro clog, you can remove and clean it or replace it with the spare suction valve secured on the mounting card.

Insertion of Endoscopic Accessories

Always make sure to select the correct size endoscopic accessory for the aScope Gastro. Inspect the endoscopic accessory before using it. If there is any irregularity in its operation or external appearance replace it. Insert the endoscopic accessory into the working channel port and advance it carefully through the working channel until it can be seen on the monitor.

- Select the endoscopic accessories compatible with the aScope Gastro. Refer to the accessories instruction manuals for operating instructions.
- Always select the correct size endoscopic accessory compatible with the aScope Gastro. There is no guarantee that instruments selected solely using this minimum working channel size will be compatible in combination. Compatibility of selected instruments should be assessed prior to the procedure.
- Confirm that the tip of the endoscopic accessory is closed or retracted into its sheath.
- Hold the endoscopic accessory approximately 4 cm from the biopsy valve and advance it slowly and straight into the biopsy valve using short strokes while observing the endoscopic image.

Withdrawal of Endoscopic Accessories

- Ensure that the accessory/tool is in a neutral position and withdraw the accessory from the aScope Gastro through the biopsy valve.
- If the accessory cannot be removed retract the aScope Gastro while observing the endoscopic image.

Withdrawal of the aScope Gastro 9

- Stop using the image magnification (zoom) function of the aBox 2.
- Aspirate accumulated air, blood, mucus or other debris by depressing the suction valve.
- Move the Up/Down angulation lock to the "F" direction to release the angulation.
- Turn the Left/Right angulation lock to the "F" direction to release the angulation.
- Carefully withdraw the aScope Gastro while observing the endoscopic image.
- Remove mouthpiece from the patient's mouth.

4.6. After Use

Detach the tubings from the Insufflation/Water, Suction and the Auxiliary Water systems from the aScope Gastro Connector. 10

Press the Disconnection Button and disconnect the aScope Gastro from the aBox 2. 11
Check the aScope Gastro for any missing parts, evidence of damage, cuts, holes, sagging, or other irregularities on the bending and insertion section including the distal end of the aScope Gastro. 12

Should any irregularities exist, immediately determine if any parts are missing and take the necessary corrective action(s).

Discard the aScope Gastro in accordance with local guidelines for medical waste with electronic components. 13

Returning Devices to Ambu

Should it be necessary to return an aScope Gastro to Ambu for evaluation, please contact your representative at Ambu beforehand for instructions and/or guidance.

To prevent infection, it is strictly forbidden to ship contaminated medical devices.

The medical device aScope Gastro must be decontaminated on site prior to shipment to Ambu. Ambu reserves the right to return contaminated medical devices to the sender.

Disposal of the aScope Gastro

The aScope Gastro is for single-use only.

Do not soak, rinse or sterilize as it may leave harmful residues or cause malfunction of the aScope Gastro. The design and material used are not compatible with conventional cleaning and sterilization procedures.

Discard all the packaging including the spare suction valve after the procedure.

5. Technical Device Specifications

5.1. Standards Applied

The aScope Gastro conforms with:

- IEC 60601-1 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance.
- IEC 60601-1-2 Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic disturbances – Requirements and tests.
- IEC 60601-2-18 Medical electrical equipment – Part 2-18: Particular requirements for the basic safety and essential performance of endoscopic equipment.
- ISO 10993-1 Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process.
- ISO 8600-1 Endoscopes – Medical endoscopes and endotherapy devices – Part 1: General requirements.

5.2. aScope Gastro Specifications

No.	Product Specification		
1	Insertion Section Dimensions		
1.1	Bending angle	Up:	210°
		Down:	90°
		Left:	100°
		Right:	100°
1.2	Max. Insertion tube diameter	10.4 mm	
1.3	Distal tip diameter	9.9 mm	
1.4	Working length	1030 mm	
2	Working Channel		
2.1	Min. working channel width	2.8 mm	
3	Optics		
3.1	Field of view	140°	
3.2	Direction of view	0° (forward pointing)	
3.3	Depth of field	3 – 100 mm	
3.4	Illumination method	LED	
4	Connections		
4.1	Medical grade CO ₂ insufflator	Max. 12 psi/80 kPa (relative pressure)	
4.2	Vacuum source	Max. -11 psi /-76 kPa (relative pressure)	
4.3	Irrigation pump	Max. 72.5 psi/500 kPa (relative pressure)	
5	Operating Environment		
5.1	Temperature	10 – 40 °C	
5.2	Relative humidity	30 – 85 %	
5.3	Atmospheric pressure	80 – 106 kPa	
6	Sterilization		
6.1	Method of sterilization	Ethylene oxide (EtO)	
7	Biocompatibility		
7.1	aScope Gastro is biocompatible		

8 Storage and Transportation Conditions		
8.1	Transportation temperature	-10 – 55 °C
8.2	Storage temperature	10 – 25 °C
8.3	Relative humidity	10 – 95 %
8.4	Atmospheric pressure	50 – 106 kPa

6. Trouble Shooting

The following table shows the possible causes of and counter measures against challenges that may occur due to equipment setting errors or damage to the aScope Gastro. Make sure to contact your local Ambu representative for detailed information.

Prior to use please do the pre-check as described in the Quick Guide.

6.1. Angulation, Bending and Angulation Locks

Possible Problem	Possible Cause	Recommended Action
Increased resistance during angulation control wheel operation.	The angulation lock is engaged.	Release the angulation lock.
Angulation lock in attracted mode.	Angulation lock is not correctly attracted.	Check whether the correct angulation lock is attracted.
Distal end does not angulate when control wheel is operated.	aScope Gastro is defective.	Withdraw the aScope Gastro and connect a new scope.
Max. bending angles cannot be achieved.	aScope Gastro is defective.	Withdraw the aScope Gastro and connect a new scope.
Bending section angulates in the opposite direction.	aScope Gastro is defective.	Withdraw the aScope Gastro and connect a new scope.

6.2. Rinsing and Insufflation

Possible Problem	Possible Cause	Recommended Action
Rinsing not possible.	Rinsing tubing not properly connected.	Connect the rinsing tubing properly to the aScope Gastro.
	CO ₂ regulator is not functioning.	Refer to the CO ₂ regulator IFU.
	Sterile water source set-up sub-optimal.	Refer to the water source IFU.
	aScope Gastro is defective.	Withdraw the aScope Gastro and connect a new scope.

Possible Problem	Possible Cause	Recommended Action
Insufflation not possible.	CO ₂ regulator is not functioning.	Refer to the CO ₂ regulator IFU.
	Sterile water source set-up sub-optimal.	Refer to the water source IFU.
	aScope Gastro is defective.	Withdraw the aScope Gastro and connect a new scope.
Continuous insufflation without operating insufflation/rinsing valve.	Insufflation/ rinsing valve opening is blocked.	Withdraw the aScope Gastro and connect a new scope.
Insufficient insufflation.	CO ₂ insufflator is not connected or not switched ON.	Connect/switch on compatible insufflator. Adjust insufflator settings.
	CO ₂ – source is empty.	Connect a new CO ₂ – source.
	Suction is activated.	Deactivate suction.

6.3. Suction

Possible Problem	Possible Cause	Recommended Action
Diminished or no suction.	Vacuum source/suction pump is not connected or not switched ON.	Connect the vacuum source/ suction pump and power ON.
	Suction canister is full or not connected.	Change the suction canister if it is full. Connect a suction canister.
	Suction valve is blocked.	Remove the valve and rinse with sterile water using a syringe and reuse the valve. Or replace with the spare suction valve.
	Biopsy valve is not properly connected.	Attach valve correctly.
	Biopsy valve cap is open.	Close cap.
	Vacuum source/suction pump to weak (min. -7 kPa).	Increase vacuum pressure.
	Vacuum source/suction pump is defective.	Replace with a new vacuum source/suction pump.
	Working channel blocked.	Flush sterile water with a syringe through the working channel.
	aScope Gastro is defective.	Withdraw the aScope Gastro and connect a new scope.

6.4. Working Channel

Possible Problem	Possible Cause	Recommended Action
Working channel or access is constricted or blocked (accessory instruments does not pass through channel smoothly).	Endoscopic accessory is not compatible.	Select a compatible accessory.
	Endoscopic accessory is open.	Close the endoscopic accessory or retract it into its sheath.
	Working channel is blocked.	Try to unblock it by flushing sterile water into the working channel with a syringe.
	Biopsy valve is not open.	Open the cap of the biopsy valve.

6.5. Image Quality and Brightness

Possible Problem	Possible Cause	Recommended Action
No video image.	aBox 2 or ancillary equipment not switched on.	Switch aBox 2 and ancillary equipment on.
	aScope Gastro connector is not properly connected to the aBox 2.	Connect the aScope Gastro connector properly to the aBox 2.
	aScope Gastro is defective.	Withdraw the aScope Gastro and connect a new scope.
	aBox 2 defective.	Contact your Ambu representative.
Image suddenly darkens.	Camera or illumination failure.	Withdraw the aScope Gastro and connect a new scope.
Blurry image.	Objective lens is dirty.	Rinse the objective lens.
	Water drops on the outside of the lens.	Insufflate and/or rinse to remove water drops from the lens.
	Condensation on the inside of the lens.	Increase the water temperature in the water container and continue to use the aScope Gastro.
	aBox 2 Image settings incorrect.	See aBox 2 IFU.
Flickering images.	Signal interference from activated high-frequency current.	Use alternative mode or settings on the HF-generator with lesser Peak Voltage (pV).

Possible Problem	Possible Cause	Recommended Action
Dark or over-illuminated image.	aBox 2 Image settings incorrect.	See aBox 2 IFU.
	aScope Gastro defective.	Withdraw the aScope Gastro and connect a new scope.
The color tone of the endoscopic image is unusual.	aScope Gastro defective.	Withdraw the aScope Gastro and connect a new scope.
Picture is frozen.	aScope Gastro defective.	Withdraw the aScope Gastro and connect a new scope.
	aBox 2 defective.	Contact your Ambu representative.

6.6. Remote Switches

Possible Problem	Possible Cause	Recommended Action
The remote switches are not working or not working properly.	aScope Gastro connector is not properly connected to the aBox 2.	Connect the aScope Gastro connector properly to the aBox 2.
	Remote switches configuration changed.	Return to standard configuration of the remote switches or change the settings.
	Wrong remote switch operated.	Operate the correct remote switch.
	aScope Gastro defective.	Withdraw the aScope Gastro and connect a new scope.
	aBox 2 defective.	Contact your Ambu representative.

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